

HRP-conjugated Affinipure Goat Anti-Rabbit IgG(H+L) antibody (HRP)

Cat#: orb688932 (CoA)

Product Name	Goat Anti-Rabbit IgG(H+L) Antibody; HRP-conjugated
Catalog Number	Orb688932
Lot Number	L0126G
Host Species	Goat
Form	Liquid
Concentration	0.1 mg/ml
Species Cross-reactivity	Rabbit
Purification method	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	ELISA, WB
Formulation	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6
Storage and Stability	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Conclusion	pass