

Avidin-Biotin-Peroxidase Complex for ELISA

Cat#: orb763209 (ELISA Manual)

SKU: orb763209

Size: 100 µl

Form Supplied: Concentrated solution (100X)

Storage Instructions: Store at -20°C for one year from date of receipt. Avoid repeated freeze-thaw cycles.

Dilution ratio: 1:100

Application

Avidin-Biotin-Peroxidase Complex (ABC-HRP) reacts with HRP substrate to produce a blue color product. The complex can be used for in vitro quantitative and qualitative analysis of enzyme-linked immunoassays.

Description

Avidin is a protein derived from both avians and amphibians that shows considerable affinity for biotin, a co-factor that plays a role in multiple eukaryotic biological processes. Avidin and other biotin-binding proteins, including Streptavidin and NeutrAvidin protein, have the ability to bind up to four biotin molecules. The bond formation between biotin and Avidin is very rapid, and once formed, is unaffected by extremes of pH, temperature, organic solvents and other denaturing agents.

Protocol

1. Bring the vial to room temperature (18-25°C) prior to use.
2. Centrifuge the vial prior to use.
3. Add 10 µl of Avidin-Biotin-Peroxidase Complex into 990 µl of PBS containing 2% BSA.
4. Add 100 µl of the prepared 1x Avidin-Biotin-Peroxidase Complex into each well. Cover with the plate sealer provided and incubate for 30 minutes at 37°C.
5. Wash the plate after the incubation.
6. Add 90 µl of HRP substrate to each well and incubate in the dark at 37°C.

Biorbyt Ltd.

5 Orwell Furlong, Cowley Road, Cambridge, Cambridgeshire
CB4 0WY, United Kingdom
Email: info@biorbyt.com | Phone: +44 (0)1223 859 353 | Fax: +44(0)1223 280 240

Biorbyt LLC

1100 Corporate Square Drive, Helix Center, Suite 221, St Louis
MO 63132, United States
Email: info@biorbyt.com | Phone: +1 (415)-906-5211 | Fax: +1 (415) 651 8558