

BIO3-BIO1 antibody

Cat#: orb784961 (CoA)

Product Code orb784961

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Uniprot No. B0F481

Immunogen Recombinant Arabidopsis thaliana BIO3-BIO1 protein Raised In Rabbit

Species Reactivity Arabidopsis thaliana

Tested Applications ELISA, WB (ensure identification of antigen)

Form Liquid

Conjugate Non-conjugated

Storage Buffer Preservative: 0.03% Proclin 300

Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Purification Method Protein A/G

Isotype IgG

Clonality Polyclonal

Product Type Polyclonal Antibody

Immunogen Species Arabidopsis thaliana (Mouse-ear cress)

Lead Time Made-to-order (14-16 weeks)

Gene Names BIO3-BIO1