


DQX1 rabbit pAb**Cat#: orb765072 (Manual)**

For research use only. Not intended for diagnostic use.


Product Name	DQX1 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Monkey
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human DQX1. AA range:551-600
Specificity	DQX1 Polyclonal Antibody detects endogenous levels of DQX1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	ATP-dependent RNA helicase DQX1
Gene Name	DQX1
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
Observed band	79kD
Human Gene ID	165545
Human Swiss-Prot Number	Q8TE96
Alternative Names	DQX1; ATP-dependent RNA helicase DQX1; DEAQ box polypeptide 1

Background similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,



Western Blot analysis of various cells using DQX1 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA).



Western blot analysis of lysates from COS7 cells, using DQX1 Antibody. The lane on the right is blocked with the synthesized peptide.