

OMG rabbit pAb

Cat#: orb765953 (Manual)


For research use only. Not intended for diagnostic use.

Product Name	OMG rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human OMG. AA range:282-331
Specificity	OMG Polyclonal Antibody detects endogenous levels of OMG protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Oligodendrocyte-myelin glycoprotein
Gene Name	OMG
Cellular localization	Cell membrane; Lipid-anchor, GPI-anchor.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal


Concentration	1 mg/ml
Observed band	48kD
Human Gene ID	4974
Human Swiss-Prot Number	P23515
Alternative Names	OMG; OMGP; Oligodendrocyte-myelin glycoprotein

Background

caution:Do not confuse oligodendrocyte-myelin glycoprotein (OMG) with myelin-oligodendrocyte glycoprotein (MOG).,function:Cell adhesion molecule contributing to the interactive process required for myelination in the central nervous system.,PTM:O-glycosylated in its Ser/Thr-rich repeat domain .,similarity:Contains 8 LRR (leucine-rich) repeats.,subunit:Binds to RTN4R.,tissue specificity:Oligodendrocytes and myelin of the central nervous system.,



Western Blot analysis of various cells using OMG Polyclonal Antibody



Immunohistochemistry analysis of paraffin-embedded human brain, using OMG Antibody. The picture on the right is blocked with the synthesized peptide.



Explore. Bioreagents.

www.biorbyt.com