


RNF113B rabbit pAb**Cat#: orb766260 (Manual)**


For research use only. Not intended for diagnostic use.

Product Name	RNF113B rabbit pAb
Host species	Rabbit
Applications	WB;IHC
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000;IHC-p 1:50-300
Immunogen	Synthesized peptide derived from the Internal region of human RNF113B.
Specificity	RNF113B Polyclonal Antibody detects endogenous levels of RNF113B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	RING finger protein 113B
Gene Name	RNF113B
Cellular localization	U2-type spliceosomal complex,
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal


Concentration	1 mg/ml
Observed band	36kD
Human Gene ID	140432
Human Swiss-Prot Number	Q8IZP6
Alternative Names	RNF113B; RNF161; ZNF183L1; RING finger protein 113B; Zinc finger protein 183-like 1
Background	similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 1 RING-type zinc finger.,



Western Blot analysis of various cells using RNF113B Polyclonal Antibody



Western Blot analysis of 293 cells using RNF113B Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).