


Galactosidase β 1L3 rabbit pAb

Cat#: orb767508 (Manual)

For research use only. Not intended for diagnostic use.

Product Name	Galactosidase β 1L3 rabbit pAb
Host species	Rabbit
Applications	IHC;IF;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human BGLL3. AA range:291-340
Specificity	Galactosidase β 1L3 Polyclonal Antibody detects endogenous levels of Galactosidase β 1L3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Beta-galactosidase-1-like protein 3
Gene Name	GLB1L3
Cellular localization	vacuole,
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
Observed band	
Human Gene ID	112937
Human Swiss-Prot Number	Q8NCI6
Alternative Names	GLB1L3; Beta-galactosidase-1-like protein 3
Background	similarity:Belongs to the glycosyl hydrolase 35 family.,



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma, using GLB1L3 Antibody. The picture on the right is blocked with the synthesized peptide.