



GPR139 rabbit pAb

Cat#: orb767589 (Manual)

For research use only. Not intended for diagnostic use.

Product Name	GPR139 rabbit pAb
Host species	Rabbit
Applications	WB;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human GPR139. AA range:181-230
Specificity	GPR139 Polyclonal Antibody detects endogenous levels of GPR139 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Probable G-protein coupled receptor 139
Gene Name	GPR139
Cellular localization	Cell membrane; Multi-pass membrane protein.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Clonality	Polyclonal



www.biorbyt.com

Concentration	1 mg/ml
Observed band	40kD
Human Gene ID	124274
Human Swiss-Prot Number	Q6DWJ6
Alternative Names	GPR139; GPRG1; PGR3; Probable G-protein coupled receptor 139; G(q)- coupled orphan receptor GPRg1; G-protein-coupled receptor PGR3
Background	This gene encodes a member of the rhodopsin family of G-protein-coupled receptors. The encoded protein is almost exclusively expressed in the central nervous system. L-tryptophan and L-phenylalanine may act as the physiologic ligands of the encoded protein. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016],



Western Blot analysis of K562 cells using GPR139 Polyclonal Antibody



Immunofluorescence analysis of LOVO cells, using GPR139 Antibody. The picture on the right is blocked with the synthesized peptide.