


PEA3 rabbit pAb

Cat#: orb768057 (Manual)

For research use only. Not intended for diagnostic use.

Product Name	PEA3 rabbit pAb
Host species	Rabbit
Applications	IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human ETV4. AA range:281-330
Specificity	PEA3 Polyclonal Antibody detects endogenous levels of PEA3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	ETS translocation variant 4
Gene Name	ETV4
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
Observed band	
Human Gene ID	2118
Human Swiss-Prot Number	P43268
Alternative Names	ETV4; E1AF; PEA3; ETS translocation variant 4; Adenovirus E1A enhancer-binding protein; E1A-F; Polyomavirus enhancer activator 3 homolog; Protein PEA3
Background	function:Transcriptional activator that binds to the enhancer of the adenovirus E1A gene; the core-binding sequence is 5'[AC]GGA[AT]GT-3'.,similarity:Belongs to the ETS family.,similarity:Contains 1 ETS DNA-binding domain.,



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using ETV4 Antibody. The picture on the right is blocked with the synthesized peptide.