


Ataxin-7L1 rabbit pAb**Cat#: orb768110 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	Ataxin-7L1 rabbit pAb
Host species	Rabbit
Applications	IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human ATXN7L1. AA range:441-490
Specificity	Ataxin-7L1 Polyclonal Antibody detects endogenous levels of Ataxin-7L1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Ataxin-7-like protein 1
Gene Name	ATXN7L1
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
Observed band	
Human Gene ID	222255
Human Swiss-Prot Number	Q9ULK2
Alternative Names	ATXN7L1; ATXN7L4; KIAA1218; Ataxin-7-like protein 1; Ataxin-7-like protein 4
Background	similarity:Contains 1 SCA7 domain.,



Immunohistochemistry analysis of paraffin-embedded human testis tissue, using ATXN7L1 Antibody. The picture on the right is blocked with the synthesized peptide.