

Olfactory receptor 5AP2 rabbit pAb

Cat#: orb768687 (Manual)


For research use only. Not intended for diagnostic use.

Product Name	Olfactory receptor 5AP2 rabbit pAb
Host species	Rabbit
Applications	WB;IF;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human OR5AP2. AA range:267-316
Specificity	Olfactory receptor 5AP2 Polyclonal Antibody detects endogenous levels of Olfactory receptor 5AP2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Olfactory receptor 5AP2
Gene Name	OR5AP2
Cellular localization	Cell membrane; Multi-pass membrane protein.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal


Concentration	1 mg/ml
Observed band	35kD
Human Gene ID	338675
Human Swiss-Prot Number	Q8NGF4
Alternative Names	OR5AP2; Olfactory receptor 5AP2

Background

Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],



Western Blot analysis of 3T3 cells using Olfactory receptor 5AP2 Polyclonal Antibody



Immunofluorescence analysis of HeLa cells, using OR5AP2 Antibody. The picture on the right is blocked with the synthesized peptide.