

## Olfactory receptor 8H3 rabbit pAb

**Cat#: orb768932 (Manual)**


For research use only. Not intended for diagnostic use.

<b>Product Name</b>	Olfactory receptor 8H3 rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	IF;ELISA
<b>Species Cross-Reactivity</b>	Human;Rat;Mouse;
<b>Recommended dilutions</b>	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human OR8H3. AA range:262-311
<b>Specificity</b>	Olfactory receptor 8H3 Polyclonal Antibody detects endogenous levels of Olfactory receptor 8H3 protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
<b>Storage</b>	Store at -20°C. Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	Olfactory receptor 8H3
<b>Gene Name</b>	OR8H3
<b>Cellular localization</b>	Cell membrane; Multi-pass membrane protein.
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal

<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	
<b>Human Gene ID</b>	390152
<b>Human Swiss-Prot Number</b>	Q8N146
<b>Alternative Names</b>	OR8H3; Olfactory receptor 8H3; Olfactory receptor OR11-172

## Background

olfactory receptor family 8 subfamily H member 3(OR8H3) Homo sapiens  
Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],



**Immunofluorescence analysis of HUVEC cells, using OR8H3 Antibody. The picture on the right is blocked with the synthesized peptide.**