


**ANKRD20 rabbit pAb****Cat#: orb769147 (Manual)**

For research use only. Not intended for diagnostic use.

|                                 |                                                                                                                       |
|---------------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Product Name</b>             | ANKRD20 rabbit pAb                                                                                                    |
| <b>Host species</b>             | Rabbit                                                                                                                |
| <b>Applications</b>             | IHC;IF;ELISA                                                                                                          |
| <b>Species Cross-Reactivity</b> | Human;Rat;Mouse;                                                                                                      |
| <b>Recommended dilutions</b>    | Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.                            |
| <b>Immunogen</b>                | The antiserum was produced against synthesized peptide derived from human ANKRD20A3. AA range:551-600                 |
| <b>Specificity</b>              | ANKRD20 Polyclonal Antibody detects endogenous levels of ANKRD20 protein.                                             |
| <b>Formulation</b>              | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..                                              |
| <b>Storage</b>                  | Store at -20°C. Avoid repeated freeze-thaw cycles.                                                                    |
| <b>Protein Name</b>             | Ankyrin repeat domain-containing protein 20A3                                                                         |
| <b>Gene Name</b>                | ANKRD20A3                                                                                                             |
| <b>Cellular localization</b>    |                                                                                                                       |
| <b>Purification</b>             | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Clonality</b>                | Polyclonal                                                                                                            |

|                                |                                                                                                             |
|--------------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Concentration</b>           | 1 mg/ml                                                                                                     |
| <b>Observed band</b>           | 95kD                                                                                                        |
| <b>Human Gene ID</b>           | 441425                                                                                                      |
| <b>Human Swiss-Prot Number</b> | Q5VUR7                                                                                                      |
| <b>Alternative Names</b>       | ANKRD20A3; Ankyrin repeat domain-containing protein 20A3                                                    |
| <b>Background</b>              | similarity:Contains 5 ANK repeats.,tissue specificity:Widely expressed; not found in kidney, nor in liver., |



**Immunohistochemistry analysis of paraffin-embedded human brain, using ANKRD20A3 Antibody. The picture on the right is blocked with the synthesized peptide.**