


pNO40 rabbit pAb**Cat#: orb769403 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	pNO40 rabbit pAb
Host species	Rabbit
Applications	IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human ZCCHC17. AA range:51-100
Specificity	pNO40 Polyclonal Antibody detects endogenous levels of pNO40 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Nucleolar protein of 40 kDa
Gene Name	ZCCHC17
Cellular localization	Nucleus, nucleolus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
Observed band	
Human Gene ID	51538
Human Swiss-Prot Number	Q9NP64
Alternative Names	ZCCHC17; PS1D; HSPC243; HSPC251; LDC4; Nucleolar protein of 40 kDa; pNO40; Pnn-interacting nucleolar protein; Putative S1 RNA-binding domain protein; PS1D protein; Zinc finger CCHC domain-containing protein 17
Background	similarity:Contains 1 CCHC-type zinc finger.,similarity:Contains 1 S1 motif domain.,subunit:Interacts with PNN. Associates with the 60S ribosomal subunit.,



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using ZCCHC17 Antibody. The picture on the right is blocked with the synthesized peptide.