

ZNF446 rabbit pAb**Cat#: orb769616 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	ZNF446 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF446. AA range:141-190
Specificity	ZNF446 Polyclonal Antibody detects endogenous levels of ZNF446 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger protein 446
Gene Name	ZNF446
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
Observed band	49kD
Human Gene ID	55663
Human Swiss-Prot Number	Q9NWS9
Alternative Names	ZNF446; ZKSCAN20; Zinc finger protein 446; Zinc finger protein with KRAB and SCAN domains 20
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 1 SCAN box domain.,similarity:Contains 3 C2H2-type zinc fingers.,