


ZNF695 rabbit pAb**Cat#: orb769744 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	ZNF695 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Immunogen	Synthesized peptide derived from the Internal region of human ZNF695.
Specificity	ZNF695 Polyclonal Antibody detects endogenous levels of ZNF695 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger protein 695
Gene Name	ZNF695
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
Observed band	65-70kD or 20kd for isoform Q8IW36-1
Human Gene ID	57116
Human Swiss-Prot Number	Q8IW36
Alternative Names	ZNF695; Zinc finger protein 695; Zinc finger protein SBZF3

Background function:May be involved in transcriptional regulation.,similarity:Contains 1 KRAB domain.,



Western Blot analysis of 293T cells using ZNF695 Polyclonal Antibody diluted at 1:2000. Secondary antibody(catalog#:RS0002) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003,Inventbiotech,MN,USA)