


Nkx-2.4 rabbit pAb**Cat#: orb770053 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	Nkx-2.4 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Monkey
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Immunogen	Synthesized peptide derived from Nkx-2.4 . at AA range: 250-330
Specificity	NKX2.4 Polyclonal Antibody detects endogenous levels of NKX2.4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Homeobox protein Nkx-2.4
Gene Name	NKX2-4
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
Observed band	36kD
Human Gene ID	644524
Human Swiss-Prot Number	Q9H2Z4
Alternative Names	NKX2-4; NKX2D; Homeobox protein Nkx-2.4; Homeobox protein NK-2 homolog D

Background function:Probable transcription factor.,similarity:Belongs to the NK-2 homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4°C over night