


CEACAM21 rabbit pAb**Cat#: orb770810 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	CEACAM21 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human CEACAM21. AA range:171-220
Specificity	CEACAM21 Polyclonal Antibody detects endogenous levels of CEACAM21 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Carcinoembryonic antigen-related cell adhesion molecule 21
Gene Name	CEACAM21
Cellular localization	Membrane ; Single-pass type I membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
Observed band	32kD
Human Gene ID	90273
Human Swiss-Prot Number	Q3KPI0
Alternative Names	CEACAM21; Carcinoembryonic antigen-related cell adhesion molecule 21

Background similarity:Belongs to the immunoglobulin superfamily. CEA family.,similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,



Western blot analysis of the lysates from HT-29 cells using CEACAM21 antibody.