

KKLC1 rabbit pAb**Cat#: orb771944 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	KKLC1 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110
Specificity	KKLC1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Kita-kyushu lung cancer antigen 1 (KK-LC-1) (Cancer/testis antigen 83) (CT83)
Gene Name	KKLC1 CXorf61
Cellular localization	Cell membrane ; Single-pass type II membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration 1 mg/ml


Observed band 12kD

Human Gene ID 203413

Human Swiss-Prot Number Q5H943

Alternative Names

Background tissue specificity:Specifically expressed in testis. Expressed by cancer cell lines.,



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4°over night