

BATF2 rabbit pAb**Cat#: orb772316 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	BATF2 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 180-260
Specificity	BATF2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Basic leucine zipper transcriptional factor ATF-like 2 (B-ATF-2) (Suppressor of AP-1 regulated by IFN) (SARI)
Gene Name	BATF2
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration 1 mg/ml


Observed band 30kD

Human Gene ID 116071

Human Swiss-Prot Number Q8N1L9

Alternative Names

Background similarity:Belongs to the bZIP family.,similarity:Contains 1 bZIP domain.,



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4°over night