

BRI3 rabbit pAb**Cat#: orb772333 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	BRI3 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	BRI3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Brain protein I3 (pRGR2)
Gene Name	BRI3
Cellular localization	Lysosome membrane ; Multi-pass membrane protein .; [Isoform 1]: Cytoplasm, perinuclear region . Co-localizes with MGAT1 and IFITM3 at the perinuclear region. .; [Isoform 2]: Cytoplasm . Nucleus . Diffuse localization in the cytoplasm and nucleus. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Clonality Polyclonal

Concentration 1 mg/ml


Observed band 13kD

Human Gene ID 25798

Human Swiss-Prot Number O95415

Alternative Names

Background similarity:Belongs to the BRI3 family.,



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night