

DBX2 rabbit pAb**Cat#: orb772592 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	DBX2 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	DBX2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Homeobox protein DBX2 (Developing brain homeobox protein 2)
Gene Name	DBX2
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration 1 mg/ml


Observed band 37kD

Human Gene ID 440097

Human Swiss-Prot Number Q6ZNG2

Alternative Names

Background similarity:Belongs to the H2.0 homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night