


GBX1 rabbit pAb**Cat#: orb772594 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	GBX1 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	GBX1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Homeobox protein GBX-1 (Gastrulation and brain-specific homeobox protein 1)
Gene Name	GBX1
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
Observed band	39kD
Human Gene ID	2636
Human Swiss-Prot Number	Q14549
Alternative Names	

Background similarity:Contains 1 homeobox DNA-binding domain.,



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4°over night