

PIGW rabbit pAb**Cat#: orb772871 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	PIGW rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 370-450
Specificity	PIGW Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Phosphatidylinositol-glycan biosynthesis class W protein (PIG-W) (EC 2.3.-.-)
Gene Name	PIGW
Cellular localization	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
----------------------	---------

Observed band	55kD
----------------------	------

Human Gene ID	284098
----------------------	--------

Human Swiss-Prot Number	Q7Z7B1
--------------------------------	--------

Alternative Names	
--------------------------	--

Background

Glycosylphosphatidylinositol (GPI) is a complex glycolipid that anchors many proteins to the cell surface. PIGW acts in the third step of GPI biosynthesis and acylates the inositol ring of phosphatidylinositol (Murakami et al., 2003 [PubMed 14517336]).[supplied by OMIM, Mar 2008],