

PKHG1 rabbit pAb**Cat#: orb772886 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	PKHG1 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 780-860
Specificity	PKHG1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Pleckstrin homology domain-containing family G member 1
Gene Name	PLEKHG1 KIAA1209
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration 1 mg/ml

Observed band 152kD

Human Gene ID 57480

Human Swiss-Prot Number Q9ULL1

Alternative Names

Background PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 DH (DBL-homology) domain.,similarity:Contains 1 PH domain.,