

DEN5A rabbit pAb**Cat#: orb772955 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	DEN5A rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 650-730
Specificity	DEN5A Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	DENN domain-containing protein 5A (Rab6-interacting protein 1) (Rab6IP1)
Gene Name	DENND5A KIAA1091 RAB6IP1
Cellular localization	Golgi apparatus membrane .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration 1 mg/ml

Observed band 141kD

Human Gene ID 23258

Human Swiss-Prot Number Q6IQ26

Alternative Names

Background

function:May be involved in Rab6-mediated GTPase signaling.,similarity:Belongs to the RAB6IP1 family.,similarity:Contains 1 dDENN domain.,similarity:Contains 1 DENN domain.,similarity:Contains 1 PLAT domain.,similarity:Contains 1 uDENN domain.,similarity:Contains 2 RUN domains.,subunit:Probably interacts with RAB6A bound to GTP.,