

TIPRL rabbit pAb**Cat#: orb773233 (Manual)**


For research use only. Not intended for diagnostic use.

Product Name	TIPRL rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TIPRL Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	TIP41-like protein (Putative MAPK-activating protein PM10) (Type 2A-interacting protein) (TIP)
Gene Name	TIPRL
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
Observed band	29kD
Human Gene ID	261726
Human Swiss-Prot Number	O75663
Alternative Names	

Background

TIPRL is an inhibitory regulator of protein phosphatase-2A (PP2A) (see PPP2CA; MIM 176915), PP4 (see PPP4C; MIM 602035), and PP6 (see PPP6C; MIM 612725) (McConnell et al., 2007 [PubMed 17384681]).[supplied by OMIM, Nov 2010],



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night