



UFD1 rabbit pAb

Cat#: orb773301 (Manual)

For research use only. Not intended for diagnostic use.

Product Name	UFD1 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	UFD1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Ubiquitin fusion degradation protein 1 homolog (UB fusion protein 1)
Gene Name	UFD1L
Cellular localization	Nucleus . Cytoplasm, cytosol .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Clonality	Polyclonal



Concentration	1 mg/ml
Observed band	33kD
Human Gene ID	7353
Human Swiss-Prot Number	Q92890
Alternative Names	
Background	The protein encoded by this gene forms a complex with two other proteins, nuclear protein localization-4 and valosin-containing protein, and this complex is necessary for the degradation of ubiquitinated proteins. In addition, this complex controls the disassembly of the mitotic spindle and the formation of a closed nuclear envelope after mitosis. Mutations in this gene have been associated with Catch 22 syndrome as well as cardiac and craniofacial defects. Alternative splicing results in multiple transcript variants encoding different isoforms. A related pseudogene has been identified on chromosome 18. [provided by RefSeq, Jun 2009],



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4°over night