

ESRP2 rabbit pAb**Cat#: orb773527 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	ESRP2 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ESRP2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Epithelial splicing regulatory protein 2 (RNA-binding motif protein 35B) (RNA-binding protein 35B)
Gene Name	ESRP2 RBM35B PP7059
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration 1 mg/ml

Observed band 79kD

Human Gene ID 80004

Human Swiss-Prot Number Q9H6T0

Alternative Names

Background epithelial splicing regulatory protein 2(ESRP2) Homo sapiens ESRP2 is an epithelial cell-type-specific splicing regulator (Warzecha et al., 2009 [PubMed 19285943]).[supplied by OMIM, Aug 2009],