

MYCL1 rabbit pAb**Cat#: orb774144 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	MYCL1 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 170-250
Specificity	MYCL1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Protein L-Myc-1 (Class E basic helix-loop-helix protein 38) (bHLHe38)
Gene Name	MYCL1 BHLHE38 LMYC MYCL
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration 1 mg/ml

Observed band 40kD


Human Gene ID 4610

Human Swiss-Prot Number P12524

Alternative Names

Background

similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein. Binds DNA as an heterodimer with MAX.,



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4^oover night