

GPR4 rabbit pAb

Cat#: orb774348 (Manual)

For research use only. Not intended for diagnostic use.

Product Name	GPR4 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 160-240
Specificity	GPR4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	G-protein coupled receptor 4 (G-protein coupled receptor 19)
Gene Name	GPR4 GPR19
Cellular localization	Cell membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration 1 mg/ml


Observed band 39kD

Human Gene ID 2828

Human Swiss-Prot Number P46093

Alternative Names

Background function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night