

MATR3 rabbit pAb**Cat#: orb774648 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	MATR3 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein AA range: 181-231
Specificity	MATR3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Matrin-3
Gene Name	MATR3 KIAA0723
Cellular localization	Nucleus matrix.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
----------------------	---------

Observed band	93kD
----------------------	------

Human Gene ID	9782
----------------------	------

Human Swiss-Prot Number	P43243
--------------------------------	--------

Alternative Names	
--------------------------	--

Background

matrin 3(MATR3) Homo sapiens This gene encodes a nuclear matrix protein, which is proposed to stabilize certain messenger RNA species. Mutations of this gene are associated with distal myopathy 2, which often includes vocal cord and pharyngeal weakness. Alternatively spliced transcript variants, including read-through transcripts composed of the upstream small nucleolar RNA host gene 4 (non-protein coding) and matrin 3 gene sequence, have been identified. Pseudogenes of this gene are located on chromosomes 1 and X. [provided by RefSeq, Aug 2013],