

LIX1 rabbit pAb**Cat#: orb774679 (Manual)**

For research use only. Not intended for diagnostic use.

Product Name	LIX1 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	LIX1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Protein limb expression 1 homolog
Gene Name	LIX1 C5orf11
Cellular localization	cytoplasm,
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration 1 mg/ml


Observed band 31kD

Human Gene ID 167410

Human Swiss-Prot Number Q8N485

Alternative Names

Background similarity:Belongs to the LIX1 family.,



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4°over night