

FAM13A antibody

Cat #: orb626914 (manual)

Product Information

Catalog No.:	orb626914
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
Isotype:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None



Immunogen information

Immunogen:	family with sequence similarity 13, member A
Synonyms:	Protein FAM13A FAM13A FAM13A1 KIAA0914
Observed MW:	65 kDa
Uniprot ID:	O94988

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200

Image

	Immunohistochemistry of paraffin-embedded human skeletal muscle using orb626914 (FAM13A antibody) at dilution of 1:50
	mouse kidney tissue was subjected to SDS PAGE followed by western blot with orb626914 (FAM13A antibody) at dilution of 1:500