


anti-Sp100 Peptide
anti-human Sp100
rabbit, polyclonal

Cat. No.	Amount
orb122509	50 µl



Description

Polyclonal Antibody pABahSp100(209-227)#2 to human Sp100.

Immunogen

The antibody binds to a peptide corresponding to amino acids 209-227 (CETEQINAKRKDTTSDKDD) of the human Sp100 protein.

Specificity

Antibody specificity was evaluated by indirect immunofluorescence and western blotting (see Figure 1A and B)

Application

Indirect Immunofluorescence, concentration: 1:1000 to 1:2000.

Western Blot, concentration: 1:500.

For *in vitro* use only

Quality guaranteed for 12 months

Store in aliquots at -20°C

Avoid freeze / thaw cycles

Form

Liquid. Supplied as serum.

BiOrbyt