

## Human IL-17RC ELISA Kit

**Catalog Number: orb138088 (Kit manual)**


Typical Data Obtained from Human IL-17RC

(TMB reaction incubate at 37°C for 20-25min)

Concentration(pg/ml)	0	250	500	1000	2000	4000	8000	16000
O.D.	0.003	0.080	0.141	0.283	0.575	1.068	1.729	2.234

Typical Human IL-17RC ELISA Kit Standard Curve

This standard curve was generated at Biorbyt for demonstration purpose only. A standard curve must be run with each assay.



- Range 250pg/ml-16000pg/ml
- Sensitivity <15pg/ml
- Specificity Natural and recombinant Human IL-17RC
- Cross-reactivity There is no detectable cross-reactivity with other relevant proteins.

### Storage

Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles(Shipped with wet ice.)

## Intra/Inter Assay Precision

**Intra-Assay Precision** (Precision within an assay) Three samples of known concentration were tested on one plate to assess intra-assay precision.

**Inter-Assay Precision** (Precision between assays) Three samples of known concentration were tested in separate assays to assess inter-assay precision.

Sample	Intra-Assay Precision			Inter-Assay Precision		
	1	2	3	1	2	3
n	16	16	16	24	24	24
Mean(ng/ml)	1.8	5	9.7	2.2	6.2	10.3
Standard deviation	0.101	0.32	0.427	0.132	0.428	0.546
CV(%)	5.6	6.4	4.4	6	6.9	5.3

## Assay Principle

Biorbyt's Human IL-17RC ELISA Kit was based on standard sandwich enzyme-linked immune-sorbent assay technology. A monoclonal antibody from mouse specific for IL-17RC has been precoated onto 96-well plates. Standards(Expression system for standard: NSO; Immunogen sequence: L21-A454) and test samples are added to the wells, a biotinylated detection polyclonal antibody from goat specific for IL-17RC is added subsequently and then followed by washing with PBS or TBS buffer. Avidin-Biotin-Peroxidase Complex was added and unbound conjugates were washed away with PBS or TBS buffer. HRP substrate TMB was used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the Human IL-17RC amount of sample captured in plate.

## Kit Components

Materials included in the kit

Description	Quantity
96-well plate precoated with anti-Human IL17RC antibody	1
lyophilized recombinant Human IL17RC standard	20ng/tubex2
biotinylated anti-Human IL17RC antibody	130ul(dilution 1:100)
Avidin-Biotin-Peroxidase Complex(ABC)	130ul(dilution 1:100)
sample diluent buffer	30ml
antibody diluent buffer	12ml
ABC diluent buffer	12ml
TMB color developing agent	10ml
TMB stop solution	10ml
Adhesive cover	4

Washing buffer Preparation: Dissolve AR0030-E to 1000ml distilled water and adjust pH to 7.2~7.6. Finally, adjust the total volume to 1L.

### Materials Required But Not Provided

1. Microplate reader in standard size.
2. Automated plate washer.

- Adjustable pipettes and pipette tips. Multichannel pipettes are recommended in the condition of large amount of samples in the detection.
- Clean tubes and Eppendorf tubes.

### Notice Before Application

Please read the following instructions before starting the experiment.

- To inspect the validity of experiment operation and the appropriateness of sample dilution proportion, pilot experiment using standards and a small number of samples is recommended.
- The TMB Color Developing agent is colorless and transparent before using, contact us freely if it is not the case.
- Before using the Kit, spin tubes and bring down all components to the bottom of tubes.
- Duplicate well assay is recommended for both standard and sample testing.
- Don't let 96-well plate dry, for dry plate will inactivate active components on plate.
- Don't reuse tips and tubes to avoid cross contamination.
- Avoid using the reagents from different batches together.
- In order to avoid marginal effect of plate incubation due to temperature difference (reaction may be stronger in the marginal wells), it is suggested that the diluted ABC and TMB solution will be pre-warmed in 37° for 30 min before using.
- Take precautionary measures to prevent operator contamination (such as saliva and other body fluids) of kit reagents while running this assay.

### Preparation

- Sample Preparation and Storage

Store samples to be assayed within 24 hours at 2-8°C. For long-term storage, aliquot and freeze samples at -20°C. Avoid repeated freeze-thaw cycles.

**Cell lysates:** After sufficient splitting, there should be no obvious cell sediment. Centrifuge cell lysates at approximately 10000 X g for 5 min. Collect the cell lysate supernates to go ahead.

**Serum:** Allow the serum to clot in a serum separator tube (about 4 hours) at room temperature. Centrifuge at approximately 1000 X g for 15 min. Analyze the serum immediately or aliquot and store samples at -20°C.

**Cell culture supernates:** Remove particulates by centrifugation, assay immediately or aliquot and store samples at  $-20^{\circ}\text{C}$ .

**Plasma:** Collect plasma using heparin or EDTA as an anticoagulant. Centrifuge for 15 min at  $1500 \times g$  within 30 min of collection. Assay immediately or aliquot and store samples at  $-20^{\circ}\text{C}$ .

## 2. Sample Dilution Guideline

The user needs to estimate the concentration of the target protein in the sample and select a proper dilution factor so that the diluted target protein concentration falls near the middle of the linear regime in the standard curve. Dilute the sample using the provided diluent buffer. The following is a guideline for sample dilution. Several trials may be necessary in practice.

**The sample must be well mixed with the diluents buffer.**

High target protein concentration ( $160000\text{pg/ml}$ - $1600000\text{pg/ml}$ ). The working dilution is 1:100. i.e. Add  $1\mu\text{l}$  sample into  $99\mu\text{l}$  sample diluent buffer.

Medium target protein concentration ( $16000\text{pg/ml}$ - $160000\text{pg/ml}$ ). The working dilution is 1:10. i.e. Add  $10\mu\text{l}$  sample into  $90\mu\text{l}$  sample diluent buffer.

Low target protein concentration ( $250\text{pg/ml}$ - $16000\text{pg/ml}$ ). The working dilution is 1:2. i.e. Add  $50\mu\text{l}$  sample to  $50\mu\text{l}$  sample diluent buffer.

Very Low target protein concentration ( $0\text{pg/ml}$ - $250\text{pg/ml}$ ). No dilution necessary, or the working dilution is 1:2.

## 3. Reagent Preparation and Storage

- A. Reconstitution of the Human IL-17RC standard: IL-17RC standard solution should be prepared no more than 2 hours prior to the experiment. Two tubes of IL-17RC standard ( $20\text{ng/tube}$ ) are included in each kit. Use one tube for each experiment.
  - a.  $20000\text{pg/ml}$  of Human IL-17RC standard solution: Add  $1\text{ml}$  sample diluent buffer into one tube, keep the tube at room temperature for 10 min and mix thoroughly.
  - b.  $16000\text{pg/ml}$  of Human IL-17RC standard solution: Add  $0.8\text{ml}$  of the above IL-17RC standard solution into  $0.2\text{ml}$  sample diluent buffer and mix thoroughly.
  - c.  $8000\text{pg/ml}$ ? $250\text{pg/ml}$  of Human IL-17RC standard solutions: Label 6 Eppendorf tubes with  $8000\text{pg/ml}$ ,  $4000\text{pg/ml}$ ,  $2000\text{pg/ml}$ ,  $1000\text{pg/ml}$ ,  $500\text{pg/ml}$ ,  $250\text{pg/ml}$  respectively. Aliquot  $0.3\text{ml}$  of the sample diluent buffer into each tube. Add  $0.3\text{ml}$  of the above

16000pg/ml IL-17RC standard solution into 1st tube and mix. Transfer 0.3ml from 1st tube to 2nd tube and mix. Transfer 0.3ml from 2nd tube to 3rd tube and mix, and so on.

Note: The standard solutions are best used within 2 hours. The 20000pg/ml standard solution should be stored at 4°C for up to 12 hours, or at -20°C for up to 48 hours. Avoid repeated freeze-thaw cycles.

- B. Preparation of biotinylated anti-Human IL-17RC antibody working solution: The solution should be prepared no more than 2 hours prior to the experiment.
  - a. The total volume should be: 0.1ml/well x (the number of wells). (Allowing 0.1-0.2 ml more than total volume)
  - b. Biotinylated anti-Human IL-17RC antibody should be diluted in 1:100 with the antibody diluent buffer and mixed thoroughly. (i.e. Add 1µl Biotinylated anti-Human IL-17RC antibody to 99µl antibody diluent buffer.)
- C. Preparation of Avidin-Biotin-Peroxidase Complex (ABC) working solution: The solution should be prepared no more than 1 hour prior to the experiment.
  - a. The total volume should be: 0.1ml/well x (the number of wells). (Allowing 0.1-0.2 ml more than total volume)
  - b. Avidin- Biotin-Peroxidase Complex (ABC) should be diluted in 1:100 with the ABC dilution buffer and mixed thoroughly. (i.e. Add 1µl ABC to 99µl ABC diluent buffer.)

## Assay Procedure

The ABC working solution TMB color developing agent and TMB stop solution must be kept warm at 37°C for 30 min before use. When diluting samples and reagents, they must be mixed completely and evenly. Standard IL-17RC detection curve should be prepared for each experiment. The user will decide sample dilution fold by crude estimation of IL-17RC amount in samples.

1. Aliquot 0.1ml per well of the 16000pg/ml, 8000pg/ml, 4000pg/ml, 2000pg/ml, 1000pg/ml, 500pg/ml, 250pg/ml Human IL-17RC standard solutions into the precoated 96-well plate. Add 0.1ml of the sample diluent buffer into the control well (Zero well). Add 0.1ml of each properly diluted sample of Human cell culture supernates, cell lysates, serum or plasma (heparin, EDTA) to each empty well. See "Sample Dilution Guideline" above for details. It is recommended that each Human IL-17RC standard solution and each sample be measured in duplicate.
2. Seal the plate with a new adhesive cover provided and incubate at 37°C for 90 min.
3. Remove the cover, discard plate content, and blot the plate onto paper towels or other absorbent material. Do NOT let the wells completely dry at any time.

4. Add 0.1ml of biotinylated anti-Human IL-17RC antibody working solution into each well, seal the plate with a new adhesive cover provided and incubate at 37°C for 60 min.
5. Wash plate 3 times with 0.01M TBS or 0.01M PBS, and each time let washing buffer stay in the wells for 1 min. Discard the washing buffer and blot the plate onto paper towels or other absorbent material. (Plate Washing Method: Discard the solution in the plate without touching the side walls. Blot the plate onto paper towels or other absorbent material. Soak each well with at least 0.3 ml PBS or TBS buffer for 1~2 minutes. Repeat this process two additional times for a total of THREE washes. Note: For automated washing, aspirate all wells and wash THREE times with PBS or TBS buffer, overfilling wells with PBS or TBS buffer. Blot the plate onto paper towels or other absorbent material.)
6. Add 0.1ml of prepared ABC working solution into each well, seal the plate with a new adhesive cover provided and incubate at 37°C for 30 min.
7. Wash plate 5 times with 0.01M TBS or 0.01M PBS, and each time let washing buffer stay in the wells for 1-2 min. Discard the washing buffer and blot the plate onto paper towels or other absorbent material. (See Step 5 for plate washing method.)
8. Add 90µl of prepared TMB color developing agent into each well, seal the plate with a new adhesive cover and incubate at 37°C in dark for 20-25min (Note: For reference only, the optimal incubation time should be determined by end user. And the shades of blue can be seen in the wells with the four most concentrated Human IL-17RC standard solutions; the other wells show no obvious color).
9. Add 0.1ml of prepared TMB stop solution into each well. The color changes into yellow immediately.
10. Read the O.D. absorbance at 450nm in a microplate reader within 30 min after adding the stop solution.

For calculation, (the relative O.D.450) = (the O.D.450 of each well) – (the O.D.450 of Zero well). The standard curve can be plotted as the relative O.D.450 of each standard solution (Y) vs. the respective concentration of the standard solution (X). The Human IL-17RC concentration of the samples can be interpolated from the standard curve.

Note: if the samples measured were diluted, multiply the dilution factor to the concentrations from interpolation to obtain the concentration before dilution.

### Summary

1. Add samples and standards and incubate the plate at 37°C for 90 min. Do not wash.
2. Add biotinylated antibodies and incubate the plate at 37°C for 60 min. Wash plate 3 times with 0.01M TBS.

3. Add ABC working solution and incubate the plate at 37°C for 30 min. Wash plate 5 times with 0.01M TBS.
4. Add TMB color developing agent and incubate the plate at 37°C in dark for 20-25min.
5. Add TMB stop solution and read.

### Background

Interleukin-17 receptor C is a protein that in humans is encoded by the IL17RC gene. By genomic sequence analysis, IL17RC gene was mapped to chromosome 3p25.3-p24.1. This gene encodes a single-pass transmembrane protein that shares limited similarity with the interleukin-17 receptor. Multiple alternatively spliced transcript variants encoding different isoforms have been detected for this gene, but the full-length nature of only three have been determined to date. The biologic activity of IL-17 is dependent on a complex composed of IL-17RA and IL-17RC, suggesting a new paradigm for understanding the interactions between the expanded family of IL-17 ligands and their receptors.

**Biorbyt Ltd.**

5 Orwell Furlong, Cowley Road, Cambridge, Cambridgeshire  
CB4 0WY, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +44 (0)1223 859 353 | Fax: +44(0)1223 280 240

**Biorbyt LLC**

Suite 103, 369 Pine Street, San Francisco  
California 94104, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +1 (415)-906-5211 | Fax: +1 (415) 651 8558