

**Human TK1 protein (C-His)**  
**Catalog Number: orb312042****Certificate of analysis**

<b>Lot #:</b>	<b>0331110</b>
<b>Expiration Date:</b>	<b>Six months from date of receipt when stored as instructed.</b>
<b>Format:</b>	<b>Liquid</b>
<b>Filtration:</b>	<b>Filtered through a 0.2µm sterile filter.</b>
<b>Storage:</b>	<b>&lt; -20°C</b>
<b>Protein Content:</b>	<b>Qualified</b>
<b>Purity:</b>	<b>99%</b>
<b>Apparent Molecular Weight:</b>	<b>28kD</b>
<b>Endotoxin:</b>	<b>&lt;1EU/ug</b>
<b>Tested by:</b>	<b>UV280</b>
<b>Tested by:</b>	<b>SDS-PAGE</b>
<b>Tested by:</b>	<b>Reducing SDS-PAGE</b>
<b>Tested by:</b>	<b>LAL</b>

**Products are for research use only. Not intended for human or animal diagnostic or therapeutic use.**