

Product Datasheet

Vazyme - ClonExpress Ultra One Step Cloning Kit (C115-02)

Description

ClonExpress Ultra One Step Cloning Kit is a new generation of recombinant cloning kits, compatible with 1 - 5 fragments homologous recombination. The unique ligase-independent system of this kit significantly reduces the self-ligation background of the vector. Restriction endonuclease sites do not need to be considered. Highly optimized 2 × ClonExpress Mix (It contains enhanced recombinase Exnase) can significantly improve the recombination efficiency and the tolerance to impurities. pCE-Zero vector is compatible with most PCR products, enabling the specific PCR products to be used directly for recombination without any treatments, which significantly simplify the procedure.

Storage

Store at -30 ~ -15°C and transport at ≤0°C.

Note

For research use only

Application notes

Simple: Ready-to-use super mix in one tube Fast: Recombination only takes 5 - 15 minutes Efficient: Efficiently clone 50 bp - 10 kb fragments, the positive rate is > 95% Compatible: Suitable for single fragment, multiple fragment and Entry cloning 'Zero' Background: Ampicillin and Kana dual resistance vector; contains CcdB gene Software support: Provide online primer design software "CE Design"

Expiration Date

12 months from date of receipt.



Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)