

Product Datasheet

Vazyme - FastPure EndoFree Plasmid Maxi Kit V2 (DC212-C1)

Catalog Number	DC212-C1
Category	Tools
Description	<p>This kit is intended for extraction from 100 - 300 ml of overnight bacterial cultures. It adopts an optimized alkaline lysis method and features a unique silica membrane adsorption technology and special Buffer ERB and Buffer ERW for specific binding of plasmids and effective removal of endotoxins, proteins, and other impurities. The kit provides high yields of plasmid DNA with minimal endotoxin residue, suitable for use in transfection of various types of cells and common experiments such as enzyme digestion, PCR, sequencing, and ligation.</p> <p>Mechanism & Workflow</p>
Usage Notes	Restricted to UK and Ireland's customers ONLY
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)