

Product Datasheet

Randox - Toltrazuril (HRP9397)

Catalog Number	HRP9397
Description	Toltrazuril - HRP Conjugate
Target	Toltrazuril
Form/Appearance	Liquid
Storage	This conjugate may be stored for up to 3 months at +2 - +8°C. For long term storage, aliquot and store at ≤-20°C.
Note	For research use only
Application notes	competitive ELISA with polyclonal antibody PAS9457/PAS9456. The conjugate is supplied as a concentrate. Only dilute the concentrate as required and use the working strength conjugate immediately after dilution.
Hazard Information	Product should be protected from light exposure. Avoid repeated freeze/thaw cycles. Avoid using azide containing buffers. Vial should be centrifuged briefly before opening to ensure all material is removed from the vial cap.
Dilution Range	1/640.
Expiration Date	12 months from date of receipt.

RANDOX
BIOSCIENCES

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)