

Product Datasheet

Vazyme - VAHTS Universal DNA Library Prep Kit for Illumina V4 (Plate) (NDB610-01)

Catalog Number	NDB610-01
Category	Tools
Description	VAHTS Universal DNA Library Prep Kit for Illumina V4 (Plate) is a 96-well plate library preparation kit for Illumina sequencers. It can convert 100 pg - 1 µg fragmented input DNA into a specialized Illumina library. The kit is applicable to prepare the PCR or PCR-Free libraries from various samples and compatible with the targeted hybridization capture. All the reagents provided in the kit have undergone rigorous quality control and functional testing to ensure the optimal stability and repeatability of library preparation.
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store at -30 ~ -15°C and ship at ≤0°C.
Note	For research use only
Expiration Date	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)