

Product Datasheet

Vazyme - VAHTS PCR-Free DNA Adapters Set 1 for MGI (NM10901)

Catalog Number	NM10901
Description	VAHTS PCR-Free DNA Adapters for MGI is a specialized kit for MGI high-throughput sequencing platform. It is suitable for preparing multi-sample PCR-Free DNA libraries for MGI high-throughput sequencing platform. It contains 96 kinds of Adapters with different indexes. All the Adapters provided in the kit have undergone rigorous quality control and functional testing to ensure the optimal stability and repeatability of library preparation.
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Apply to MGI platform for PCR-Free DNA library preparation. Providing 96 different single indexes.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)