

Product Datasheet

Vazyme - dATP (100 mM) (P034-01)

Catalog Number	P034-01
Category	Tools
Description	dUTP/dATP/dTTP/dCTP/dGTP are high-purity single component products. The concentration of products was 100 mM (100 mmol/L). It is applicable for PCR, real-time PCR, RT-PCR, cDNA or common DNA synthesis, primer extension reaction, DNA sequencing, DNA labeling and other conventional molecular biology reactions.
Application notes	PCR companion: Individual components for selection, configurable with different concentrations of dNTPs Quality control: >99% purity confirmed by HPLC; free from qPCR, PCR, and reverse transcription inhibitors
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)