

Product Datasheet

Vazyme - AceQ Universal U+ Probe Master Mix V2 (Q513-03)

Catalog Number	Q513-03
Category	Tools
Description	AceQ Probe Master Mix is a special reagent for qPCR using probe method. The core component AceTaq DNA Polymerase is a chemically modified hot-start DNA polymerase. Combined with optimized buffer for qPCR, it can effectively inhibit non-specific amplification, thereby significantly improving the amplification efficiency. It is suitable for high-sensitivity qPCR. This product is a 2 × master mix. It can obtain good standard curves in a wide quantitative area, and the target gene can be accurately quantified and detected, with good repeatability and high reliability.
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store at -30 ~ -15°C and transport at ≤0°C. Keep away from light.
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)