

Product Datasheet

Vazyme - ChamQ Geno-SNP Probe Master Mix (Q811-03)

Catalog Number Q811-03

Description ChamQ Geno-SNP Probe Master Mix is specially designed for single nucleotide polymorphism (SNP) typing by probe-based qPCR, which can be performed simply by only adding additional primers, Taqman MGB probes and templates. This master mix uses Champagne Taq DNA Polymerase as the core enzyme, with carefully optimized Buffer, the success rate of typing on low-concentration templates and complex templates has been increased. This product contains a unique ROX Passive Reference Dye that is suitable for all qPCR instruments. The concentration of ROX does not need to be adjusted on different instruments

Storage Store at -30 ~ -15°C and transport at ≤0°C. Keep away from light.

Note For research use only

Application notes Compatible with direct amplification of blood lysate: SNP typing can be performed using only blood lysate as the template. Broad template input compatibility: It is compatible with template input ranging from 1 - 10 ng. Great GC-rich template compatibility: It accurately genotypes at SNP sites with a GC content of 25 - 73%.

Expiration Date 12 months from date of receipt.



Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)