

Product Datasheet

Vazyme - Animal Detection U+ Probe qPCR Super PreMix (QV114-02)

Catalog Number	QV114-02
Category	Tools
Description	Animal Detection Probe qPCR Super PreMix is a specialized master mix for probe qPCR to detect African swine fever virus (ASFV). Only need to add extra primers, probes and templates into the reaction system. This master mix employs an upgraded hot start Taq DNA polymerase with a carefully optimized buffer to increase detection sensitivity for low-concentration templates, while supporting the primer probe to be premixed in advance. This product can eliminate the contamination of DNA product carryover, ensuring the accuracy of the results. It is also compatible with fast program, shortening test time.
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)