

Product Datasheet

Vazyme - HiScript III All-in-one RT SuperMix Perfect for qPCR (R333-00)

Catalog Number	R333-00
Category	Tools
Description	<p>HiScript III All-in-one RT SuperMix Perfect for qPCR is an upgraded version of HiScript III RT SuperMix for qPCR (+gDNA wiper). Within one step, genomic removal and reverse transcription reactions can be achieved simultaneously. The simpler and faster operation ensure reduce risk of sample contamination and RNA degradation. This product contains HiScript III Reverse Transcriptase, heat-labile DNase and an optimized buffer system. The heat-labile DNase works quickly and efficiently; and it could be inactivated easily. The RT product is compatible with dye-based and probe-based qPCR, and can perform high-performance gene expression analysis.</p>
Usage Notes	Restricted to UK and Ireland's customers ONLY
Storage	Store at -30 ~ -15 °C and transport at ≤0 °C .
Background	Only available for the UK and Republic of Ireland
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)