

Product Datasheet

Human CD4 Mouse Monoclonal Antibody (APC) (orb1000116)

Catalog Number	orb1000116
Description	Human CD4 Mouse Monoclonal Antibody (APC)
Species/Host	Mouse
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human CD4
Target	CD4
Preservatives	0.01M PBS (pH7.4), 0.5% rAlbumin, 0.02% Proclin300.
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
Clone Number	HIT4a
Antibody Type	Primary Antibody
Dilution Range	Flow-Cyt=10ul/Test
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)