

Product Datasheet

Phospho-Mcl1 (Ser159 + Thr163) Rabbit Polyclonal Antibody (Cy7) (orb1047989)

Catalog Number	orb1047989
Category	Antibodies
Description	Phospho-Mcl1 (Ser159 + Thr163) Rabbit Polyclonal Antibody (Cy7) is a Cy7 conjugated antibody targeting MCL1. This antibody is suitable for FC, ICC, IF. It exhibits reactivity with Human, Rat samples.
Target	MCL1
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Cy7
Reactivity	Human, Rat
Predicted Reactivity	Bovine, Canine, Equine, Mouse
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Mcl 1 around the phosphorylation site of Ser159/Thr163 DG(p-S)LPS(p-T)PP
Tested applications	FC, ICC, IF

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Dilution range	Flow-Cyt=1ug/Test, IF=1:100-500, ICC/IF=1:100-500
Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)