

Product Datasheet

Phospho-Glucocorticoid Receptor (Ser211) Rabbit Polyclonal Antibody (Cy7) (orb1053160)

Catalog Number	orb1053160
Description	Phospho-Glucocorticoid Receptor (Ser211) Rabbit Polyclonal Antibody (Cy7)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Cy7
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Glucocorticoid Receptor around the phosphorylation site of Ser211 NE(p-S)PW
Target	NR3C1
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Dilution Range	Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)